

# WARNING

## Potential for Incorrect Blood Glucose Reading

Patient Name

Is using **EXTRANEAL** (Icodextrin)  
peritoneal dialysis solution

**ONLY** use glucose monitors and test strips that are glucose-specific. These methods are common in clinical laboratories. Contact the manufacturer of the glucose monitors and test strips to determine the method that is used.

**DO NOT** use glucose monitors or test strips that utilize glucose dehydrogenase pyrroloquinolinequinone (GDH-PQQ) or glucose-dye-oxidoreductase methods. In addition, some but not all monitors or test strips that utilize a glucose dehydrogenase flavin-dinucleotide (GDH-FAD) method should not be used. Use of these methods may result in falsely elevated blood glucose readings in patents using **EXTRANEAL** (icodextrin) due to maltose interference. Falsely elevated blood glucose readings may mask true hypoglycemia or lead to the erroneous diagnosis of hyperglycemia, leading to life-threatening events. Visit [www.glucosesafety.com](http://www.glucosesafety.com) for additional information, including a glucose monitor compatibility list.

**Baxter**

PSTGL03UK3 Date of Preparation May 2010. ADV10/1546R. **BAXTER** and **EXTRANEAL** are trademarks of Baxter International Inc.