

Ångelholm  
April 10, 2026\*

HemoCue AB  
Kuvettgatan 1  
SE-262 71 Ångelholm  
Sweden

To users of HemoCue® Glucose 201 Microcuvettes (Vial)

### Field Notice

**Identifier: FA-0051**

**Affected product: HemoCue® Glucose 201 Microcuvettes (Vial)**

**Article number: 110706, 110716, 110719, 110720, 110727**

**Affected LOT numbers:**

2601239  
2601244  
2601256  
2601257  
2602279  
2602294  
2602306  
2601259  
2602307  
2602308  
2601255  
2602272  
2602292  
2602273  
2603323  
2603332  
2603324  
2603325

**Unique Device Identifier (UDI-DI):**


07311091107060, 07311091107169, 07311091107190, 07311091107206,  
07311091107275

**Dear Customer**

HemoCue is issuing this Field Notice to inform customers of a potential lid orientation deviation affecting a limited number of HemoCue® Glucose 201 Microcuvette vials. This notification is provided for awareness only and does not require product return. This notification is being issued to ensure correct handling of the product with this issue and to minimize the risk of improper storage after opening.



HemoCue assures you that product integrity is not compromised when the lid is opened and closed completely. If you receive a vial where the lid appears incorrect (see picture below), please note that when opening the vial the whole lid may detach. If this happens, please take care to reattach the lid to the vial and ensure it is closed and resealed properly. As long as the lid is placed back properly, the product quality will remain maintained.




Our aim is to ensure transparency and continued safe and effective use of our products.

#### **Transmission of the Field Notice**

Our records indicate you have received one or more of the affected batches. This notice needs to be passed on to everyone within your organization that needs to be aware of this information. HemoCue kindly requests your cooperation in completing and returning the enclosed Field Notice Verification Form to confirm your receipt of this letter, regardless of whether you have product in your inventory or not.

We apologize for any inconvenience this may cause. If you have any questions regarding this Field Notice, please send an e-mail to [cuvettes@hemocue.se](mailto:cuvettes@hemocue.se).

Yours sincerely,  
HEMOCUE AB



Manuel Gonzalez Artunduaga  
Global Product Manager, Diabetes



Brian Bratz-Boyd  
Senior Director RA/QA

\*Typographical error of April 10, 2025 has been corrected to April 10, 2026 on April 17<sup>th</sup>, 2026

## Field Notice Verification form

**Identifier:** FA-0051

I hereby confirm that I have read and understood the information in this Field Notice concerning HemoCue® Glucose 201 Microcuvettes (Vial).

I confirm that this notice has been passed on to all that need to be aware within my organisation.

<b>Date:</b>	
<b>Your Name:</b>	
<b>Signature:</b>	
<b>Email Address:</b>	
<b>Customer name:</b>	
<b>Full address:</b>	

**Please sign and return, no later than May 15<sup>th</sup>, 2026, to your local distributor.**

